

LiquiCapture™ MT



An easy, cost-effective solution to collect urine samples, incorporate on-site screening, and send samples to the lab for confirmatory testing. Ensure accuracy and minimize risk of sample contamination.

- Collection Cup easily integrates 5-panel Drug Test Dip Cards for on-site screening*

**Test dip cards not included.*

Easy and quick to use.



1

Flip cup upside down and return to upside position.



2

Open flip top lid, ensure chamber is liquid filled.



3

Place testing card in open slot.



4

Confirmation needed? Peel off sticker to vacutainer chamber.



5

Depress vacuum tube into chamber to collect liquid



6

Label and store collected sample for bulks to lab.



• 1043 Pedigo Way, STE 15 Bowling Green KY 42103
• Phone: 270-904-5097 • Fax: 270-904-5864
• acculabtox.com